

VIRUS ARMOR

UV-C Disinfection Bag

www.virusarmor.com

4924 Cambridge Rd,
Fort Worth, Texas, 76155, USA
1 855-929-8879

orders@virusarmor.com



UV-C DISINFECTION BAG
EPA REGISTERED

UV-C DISINFECTION BAG



Features

- Kills Viruses, Bacteria & Germs.
- UV-C 270+/-nm
- Sterilizes in 2 minutes.
- Built in organizer.

Available in **Pink** & **Gray**



Side storage
for charging cord



Convenient & Portable

Protects you from Viruses, Bacteria & Germs Anytime, Anywhere!



Specifications

- High disinfection rate.
- Sterilizes in 2-5 minutes.
- Durable, UV LED lights.
- Lifespan 10,000 hours.
- Quality Quartz LED emits high-intensity UV rays.
- Reflective material inside enhances sterilization effect.
- Safety protection, Automatic shut off when bag is opened.
- Built-in rechargeable Lithium battery.
- Easy Disinfection, One-button Start.
- Convenient & Portable.

WARNING & SAFETY

DANGER: Risk Group 3 Ultra Violet product. These lamps emit high power UV radiation that can cause severe injury to skin and eyes. Avoid eye and skin exposure to unshielded product. Use only in an enclosed environment that shields users from the radiation.